


VF-800-0Z

Technical specs



PRECISIONARY



| | |
|----------------------------|---|
| Advance Speed | 0-20 mm / s, adjustable |
| Return Speed | 20 mm / s, fixed |
| Vibration Frequency | 0-20 Hz, adjustable |
| Vibration Amplitude | 0.4 mm, fixed |
| Z-axis Vibration | < 4 μ m |
| Blade | tungsten carbide |
| Cutting Angle | 20 degrees, fixed (standard) |
| Thickness Adjustment | digital/automatic |
| Slices Thickness | Minimum: 40 μ m Maximum: 500 μ m |
| Slice Thickness Resolution | 2 μ m |
| Tissue Diameters | 75mm, 100mm, 125mm, 150mm |
| Maximum Tissue Length | 150 mm |
| Cutting Mode | single/continuous selectable |
| Bath | 350 x 180 x 45 mm |
| Power Source | AC 100-110 V, 200-240 V |
| Power Consumption | 8W |
| Dimension | 480 x 480 x 450 mm |
| Weight | 15 kg |

Precisionary Instruments LLC
90 Canal Street, Suite 400
Boston, MA 02114

Tel: 617-682-0586
Fax: 1(866)424-2217
Email: info@precisionary.com