

VF-200

Technical specs



PRECISIONARY



| | |
|-------------------------|--|
| Advance Speed | 1-20 mm / s, adjustable |
| Return Speed | 20 mm / s, fixed |
| Vibration Frequency | 0-20 Hz, adjustable |
| Vibration Amplitude | 2 mm, fixed |
| Z-axis Vibration | < 2 μ m |
| Blade | double edge razor blade |
| Cutting Angle | 13 degrees, fixed (standard) |
| Thickness Adjustment | manual |
| Micrometer Resolution | 10 μ m / div |
| Maximum Tissue Diameter | 12.5 mm (standard tube), 15.5 mm (large tube), 6 mm (small tube) |
| Maximum Tissue Length | 25 mm |
| Minimum Slice Thickness | 10 μ m |
| Cutting Mode | Single |
| Bath | 140 x 60 x 30 mm |
| Power Source | DC 13-15 V |
| Power Consumption | 500 mW |
| Dimension (L x W x H) | 360 x 210 x 190 mm |
| Weight | 5 kg |

Precisionary Instruments LLC
90 Canal Street, Suite 400
Boston, MA 02114

Tel: 617-682-0586
Fax: 1(866)424-2217
Email: info@precisionary.com