




PRECISIONARY



Acute cortical slice for electrophysiology from 18-month old mouse brain sliced with the Compressome®



Immuno-stained astrocytes from fixed mouse brain sectioned with the Compressome®

Introducing the **NEW** vibrating microtome **Compressome® VF-310-0Z**

- 01 **EASY**
to Use
- 02 **EASY**
to Clean
- 03 **MORE**
Control Features
- 04 **BETTER**
Durability

Durable Vibrating Head Cover

Slice Number & Total Thickness Counters

- Track # of slices cut
- Displays total thickness (µm) cut

Blade Holder


- No need to align blades
- Easy to attach/remove

Auto Zero-Z® Technology

- No calibration needed
- Eliminates z-axis movement



Slice Thickness
Thickness range: 4µm - 1000µm




Buffer Tray
Easily removable

Tube Sizes
Fits up to 20mm diameter

Compression Technology
Produces even, consistent tissue slices

Buffer Tray Valve
Easy open/close to drain



Contact Us Today for a **FREE** Instant Quote!